CERAMCO® PFZ PRESS

Directions For Use
Mode d'emploi
Istruzioni per l'uso
Gebrauchsanleitung
Instrucciones de Uso

CERAMCO® PFZ PRESS

PORCELAIN FOR ZIRCONIA

DENTSPLY
CERAMCO

DENTSPLY International, Inc.
570 West College Ave.
York, PA 17405-0872
717-699-4190
www.ceramco.com
www.dentsply.com


EC REP DeguDent GmbH
Rodenbacher Chaussee 4
63457 Hanau-Wolfgang
GERMANY
+49/6181/59-50

CE R^{only}
0086

DENTSPLY
CERAMCO

The chemistry of a great smile.™


Figure 1

Figure 2

Figure 3

Figure 10


Figure 11

Figure 12

Figure 4

Figure 5

Figure 6

Figure 13

Figure 14

Figure 15

Figure 7

Figure 8

Figure 9


Figure 16

Figure 17

Figure 18

CERAMCO® PFZ PRESS

PORCELAIN FOR ZIRCONIA

PORCELAIN FOR ZIRCONIA INSTRUCTIONS

Pressable Technique

*(For Porcelain Build-Up Technique
Directions for Use consult REF 090862.)*

The Ceramco PFZ Press ingots are designed to be used only with Cercon and other 100% dense zirconia frameworks with a CTE of ~10.5 ppm/°C. Ceramco PFZ Press Paste Liner must be used with Ceramco PFZ Press Dentin Ingots. Ceramco PFZ should be used as a veneering porcelain for cut-back technique.

Indications

The Ceramco PFZ Press may be used for all single and multiple units Cercon and other 100% dense zirconia frameworks with a CTE of ~10.5 ppm/°C.

Contraindications

The Ceramco PFZ Press is not indicated for any other all-ceramic cores other than zirconia as indicated above. The Ceramco PFZ Press is also not indicated for any alloy frameworks.

Warnings/Adverse Reactions

Consult Ceramco PFZ MSDS (Material Safety Data Sheets) for safe handling and usage.

Precautions

Ceramco PFZ Powder Liner should not be used with the Ceramco PFZ Pressable Technique. Refer to the special notes throughout the instructions for best practices for using Ceramco PFZ materials.

• Porcelain Furnace

For optimum results, be sure your porcelain furnace is accurately calibrated and manufacturer recommended firing times and temperatures are followed. If necessary, adjust the furnace parameters to obtain the correct visual indicator. Ensure the furnace muffle is not contaminated with silver. Ceramco PFZ can discolor like other porcelains.

• Sagger Trays

During the firing cycle, large sagger trays will absorb some of the heat and may affect the appearance of the porcelain. For optimum results, only use the honeycomb style trays with the all-ceramic metal firing pins.

• Liquids

Ceramco PFZ liquids are recommended for their ability to improve the handling and carving of porcelain and to reduce slumping during build-up. Distilled water may also be used.

• Storage

No special storage is required. The lot number is printed on each package; please refer to the lot number in all correspondence.

Symbols Key:

Consult MSDS

Harmful

For professional dental use only

Flammability

Imperial/Metric Conversion Key:

1.0 oz. = 28.4 g	0.07 fl. oz. = 2 ml
0.1 oz. = 3.0 g	0.50 fl. oz. = 15 ml
	3.50 fl. oz. = 100 ml

Ceramco® PFZ Margin Porcelain

Note: The porcelain margins are applied and fired prior to the liner application. For margin porcelain application, see instructions beginning on page 5.

Ceramco® PFZ Press Paste Liner Application

Note: Only Ceramco PFZ Press Paste Liner can be used with the Ceramco PFZ Press Dentin Ingots. Do not use Ceramco PFZ Powder Liner with Ceramco PFZ Press Dentin Ingots.

1st Paste Liner Application

1. Ensure the zirconia framework is clean and dry. If necessary, blast the zirconia framework with aluminum oxide and then steam clean.
2. Select the proper liner depending on the shade. Express a small amount of paste from the syringe and apply with a flat brush to the coping.
3. Achieve an even coat and do not allow the material to puddle. Excessive thickness of the liner layer will fissure during firing. Lightly condense and dry the liner using a warm muffle. Always brush any dried liner that may have accumulated inside the restoration away at this time.
4. Dry and fire according to the recommended temperatures. The correct visual indicator is a slight shiny surface.

Figure 1. Clean substructure

Figure 2. First liner application

2nd Paste Liner Application

If necessary apply a second layer of paste with a flat brush, dry and fire the restoration according to the recommended firing temperatures. The correct visual indicator is a slightly shiny surface.

Wax-Up Procedures

Wax the restoration using a conventional inlay wax.

Note: Use a non-contaminating inlay wax. After burning out, some inlay waxes may leave an ash or carbon residue. This will result in dark areas on the pressed crown. Recreate all necessary anatomical features in wax. Be sure the wax has a minimum thickness of 0.8 mm. If the coping is 0.5mm, the normal total thickness should be at least 1.5 mm.

Figure 3. Wax to full contour

Figure 4. Minimum wax thickness of 0.8mm.

Spruing Procedures Determining How Many Ingots to Use

1. Weigh the sprue base without the wax patterns.
2. Attach the wax patterns to the sprue base and reweigh the sprue base.
3. The difference between the two is the weight of the wax.
4. If the weight of the wax is 0.0 - 0.6gm or less, only one 2 gram ingot is required.
5. If the weight of the wax is 0.61 - 1.4gm or less, two 2 gram ingots are required.

Attaching the Sprues

1. Select either a 100g or 200g sprue base. Be sure the sprue base is clean and free of any debris.
2. All sprues should be 2-3 mm in diameter or 8 gauge round. Extremely large cases may require a 6 gauge round sprue.
3. Sprues should be 3-5 mm in length. Do not use reservoirs.
4. A single sprue is normally attached directly to the incisal edge or on a cusp tip.
5. Do not flare the wax at the point of attachment or use sharp angles.
6. Attach the sprues to the sprue base in the same manner.
7. Allow a minimum of 3.0 mm between each wax-up, and between the wax-ups and the exterior wall of the investment.
8. Allow a minimum of 8.0 mm between the top of the wax-up and the ring positioner.

Figure 5. Proper sprue length

Figure 6. Crown place on the sprue base

Investing Procedures

Mixing the Investment

1. The 100gm ring will require one package of investment. The 200gm ring will require 2 packages of investment. Ceramco Press Investment Powder(530301) and Ceramco Press Investment Liquid(530302) are recommended.
2. For each 100gm packet of investment use 20 ml of investment liquid and 5 ml of distilled water.
3. Place the liquid in a dry mixing bowl followed by the investment.
4. Hand-mix the materials for 15 seconds.
5. Use a full vacuum and mix the materials under slow speed for 60 seconds.
6. Carefully paint the restorations with the investment and pour the remaining investment into the ring.
7. Place the ring height positioner on the top of the ring. Excess investment will push through the top of the positioner. This step is important, as it will determine the proper height of the investment for the pressing furnace.

Figure 7. 100g Sprue Set
Figure 8. Paint the crown and fill the ring

Burnout Procedures

Preparing the Invested Ring for Burnout

1. Allow the investment to bench set a minimum of 20 minutes.
2. Remove the ring height positioner and the sprue former from the investment.
3. Carefully remove any rough investment from the side and the top of the ring. Do not decrease the height of the ring. Be sure to maintain the 90° angle between the ends and the sides of the ring.
4. Do not allow any investment or debris to enter the sprue hole. If the ring is allowed to set over one hour, place the ring in a humidior until ready to burnout.

Option #1: Used when the investment has bench set only 20 min

1. Preheat the burnout furnace to 850°C
2. Place the rings directly into the preheated furnace.
3. Be sure the floor of the burnout furnace has ridges to allow the wax to easily escape. Placing the ring in the burnout furnace on a flat surface may cause the ring to crack.
4. Place the alumina plungers into the preheated furnace.
5. Allow the rings to heat soak 45 min.
6. Add 15 minutes for each additional ring.

Option #2: Used when the investment has bench set over 30 minutes

1. Place the rings and alumina plungers into a burnout furnace at room temperature.
2. Be sure the floor of the burnout furnace has ridges to allow the wax to easily escape. Placing the ring in the burnout furnace on a flat surface may cause the ring to crack.
3. Raise the temperature of the furnace at 15°C per min to 850°C.
4. Allow the rings to heat soak 1 hr.
5. Add 10 minutes for each additional ring.

Figure 9. Burning out

Pressing Procedures

Use program as indicated on Pressing chart page 9. Preheat the plungers only. Do not preheat the ingots.

Figure 10. Insert the ingots and the plunger
Figure 11. Completed pressing

Divesting and Finishing Procedures

Divesting Procedures

1. Take an unused alumina plunger and place along side the pressed ring.
2. Use a pencil and mark the approximate height of the used alumina plunger.
3. Use a disk or 50 micron aluminum oxide and cut the investment down to the top of the plunger. Carefully break off the used plunger. The plunger should be thoroughly cleaned with aluminum oxide and used again.
4. Use 50 micron aluminum oxide with 40 psi to remove the bulk of the investment. Do not blast the pressed ceramic with the aluminum oxide. The aluminum oxide is too abrasive and will harm the pressed ceramic.
5. Use 50-micron glass beads with 20 psi to remove the investment around the pressed glass. Be sure to remove all the investment.

Figure 12. Divested crown

Finishing Procedures

1. Conventional abrasives used for ceramics may be used for the Ceramco PFZ Press. Use a thin diamond disk to remove the sprues from the ceramic crowns.
2. Use a diamond bur or suitable aluminum oxide stones to recontour the sprue attachment. Do not use abrasives, which create excess heat. This will cause a thermal crack and the ceramic coping will have to be remade.
3. Cutback the incisal area to make room for the enamel porcelain. Do not grind through the dentin and expose the liner and/or zirconia coping.

Figures 13-15. Cutting back for the enamel porcelain

Porcelain Application

Note: When firing Ceramco PFZ Press crowns with PFZ glaze or PFZ porcelain it is important to use a honeycomb sagger tray with all-ceramic firing pins. These pins are easy to bend and do not act as a heat sink. Other firing pins will not allow the proper heat transfer and will create a stress point.

1. Use 50-micron aluminum oxide at 20 psi and

lightly blast the exterior surface of the ceramic coping. Be careful not to harm the margins.

2. Use a steam cleaner to clean the surfaces or place the restorations in distilled water and clean in an ultrasonic cleaner for 10 minutes.
3. Apply and fire the Ceramco PFZ dentin and/or enamel to complete the restorations. All the components of the Ceramco PFZ system may be used to enhance the esthetics and contours of the restoration.
4. Refer to the Ceramco PFZ instructions for specific details and firing instructions.
5. When the Ceramco PFZ porcelain and the Ceramco PFZ Press pressed ceramic are both on the exterior surface, the stain and/or glaze is applied in the normal manner. However after firing the pressed ceramic surface will have a slightly rough and textured surface. This surface may either be mechanically polished, or a second layer of stain and/or glaze may be applied and fired.

Figure 16. Porcelain application

Figure 17. Stain and glaze

Figure 18. Shade verification

CERAMCO® PFZ

PORCELAIN FOR ZIRCONIA

PORCELAIN APPLICATION INSTRUCTIONS

Ceramco® PFZ Margin Porcelain

Note: *The porcelain margins are applied and fired prior to the liner application.*

Tooth Margin Preparation

A 90° shoulder preparation or deep chamfer is recommended. Normal tooth reduction (1.5mm) is recommended.

Substructure Design

Extend the facial surface of the zirconia substructure so that the zirconia ends at the axial gingival corner of the margin preparation or extends slightly onto the gingival ledge.

Die Sealer Application

Using the Ceramco PFZ Die Sealer, apply a thin, even layer of sealer over the die margin area. Blow off any excess die sealer and allow the liquid to thoroughly dry. Reapply the die sealer at least 2 more times and allow to dry. Mark the facial margin with a wax pencil.

Die Release Application

Using the Ceramco PFZ Die Release, apply a thin even layer of die release over the margin area. Blow off the excess die release. Apply and allow to dry several more times until the die has a very slight sheen.

First Margin Application

Note: *Dispense only enough material for immediate use to avoid waste. Dried Ceramco PFZ margin porcelain can not be reused. Select the correct margin porcelain from Table 1.*

1. Place a small amount of the shaded margin porcelain on a clean palette. If modifier is necessary, add the modifier porcelain to the shaded margin porcelain.
2. Mix the Ceramco PFZ Margin Liquid with the margin porcelain to a thick creamy consistency. This mixture will air harden in approximately 5 minutes. If additional working time is required, mix in a small amount of distilled water.

3. Place the porcelain mixture in the gingival area, pushing the mixture down to but not over the gingival margin. Do not over build.

Note: *All porcelains shrink by volume. The larger the porcelain mass the more apparent the firing shrinkage. Better results are achieved by applying a small amount of porcelain on the first porcelain application.*

4. Allow the Ceramco PFZ margin porcelain to dry completely (at least 5 minutes). The porcelain surface will be chalky and hard. Carefully lift the coping from the die.
5. Dry and fire the restoration according to the recommended temperatures. The correct visual indicator will be shiny yet grainy.
6. Complete the 3 minute dry and 3 minute preheat cycles to remove organic material and avoid discoloration. For furnaces without a "preheat", dry the porcelain for 6 minutes. Fire the porcelain on honeycomb sagger trays with all-ceramic metal firing pins using the recommended firing parameters.

Second Margin Application

1. Be sure the die is clean and reapply the die release. Allow to dry.
2. Place the coping on the die. Make a new mix of margin porcelain and Ceramco PFZ Margin Liquid. Apply the mixture to the gingival margin.
3. Allow the margin porcelain to dry completely (at least 5 minutes) and remove the coping from the die.
4. Dry and fire the restoration according to the recommended temperatures. The correct visual indicator will be shiny yet grainy.

Note: *If additional applications are required, repeat the procedures used for the second margin application.*

5. Apply and fire the liner according to the recommended procedures.

Ceramco® PFZ Porcelain Application

Ceramco PFZ Opaqueous Dentins

Natural teeth exhibit numerous color zones. These areas have various degrees of translucency, chroma and opacity. The experienced technician learns to use a variety of shades and modifiers to effectively mimic these characteristics. Frequently, the task is compromised by circumstances that do not allow for adequate porcelain thickness. The Ceramco PFZ Opaqueous Dentin porcelains address the problem of limited space, simplify internal color modifications and aid the overall appearance of the shade.

Applications

The Ceramco PFZ Opaqueous Dentin porcelains are approximately the same hue as the Ceramco PFZ dentin shades, but have approximately 10% more opacity. They may be used undiluted, or mixed with the Ceramco PFZ dentins or dentin modifiers.

Thin Areas

Frequently insufficient dentin/incisal thickness exists (less than 1.0mm) over the liner porcelain. The shade in this thin area can be improved by applying a layer of the opaqueous dentin porcelain under the same dentin porcelain shade. When contouring, you will be able to reduce the area where the opaqueous dentin was layered and the prescribed shade will be maintained. Complete the building and firing of the restoration in the normal manner.

Pontics

A color imbalance frequently occurs between the pontics and the abutments of a multiple-unit fixed partial denture. The color difference is caused by the absence of a full zirconia substructure in the pontic area. Ceramco PFZ opaqueous dentin porcelain applied to the bottom and the cervical area of the pontic creates a uniform color completely across the fixed partial denture.

Lingual Fossas and Occlusal Surfaces

Often, inadequate space is available on the lingual surface of an anterior unit, or on the occlusal surface of a posterior unit. Ceramco PFZ opaqueous dentin porcelain applied alone or as a blend with 10% to 20% dentin modifiers in these thin areas prevents the opaque porcelain from being highlighted.

Ceramco PFZ Dentin and Dentin Modifier Porcelain

The Ceramco PFZ Dentin Modifiers have the same hue and approximately 40% less opacity than the Ceramco PFZ dentin porcelain but are approximately 70% higher in chroma (color saturation). The dentin modifiers can be used undiluted or mixed with the Ceramco PFZ Dentins, Opaqueous Dentins, Natural Enamels, or Opal Enamels for a variety of effects. Select the recommended dentin modifier that corresponds with the desired shade according to the following chart. For a very subtle effect, mix in a very small amount of dentin modifier with the dentin porcelain. Select the correct dentin modifier using **Table 1**.

Ceramco PFZ Natural Enamels, Opal Enamels, and Mamelon Porcelains

The Ceramco PFZ Natural Enamels are non-opalescent and can be used in thin or thick layers over the dentin to duplicate the vitality of natural teeth. They can also be mixed with the dentin or opal enamel porcelains for a variety of effects. The Ceramco PFZ Opal Enamels are an intense opal enamel porcelain. They can be mixed in with the Natural Enamels or used as an effect powder. The mamelon porcelains can be used in full concentration or diluted with the clear enamel porcelain. They are normally placed on the dentin porcelain lobes prior to the enamel porcelain application.

Basic Porcelain Layering Procedures

1. Areas on the model that the porcelain will contact should be completely sealed with Ceramco PFZ Die Sealer. Other sealers or oils may leave a residue that will cause discoloration and should not be used.
2. Mix the porcelains with Ceramco PFZ Modeling Liquid, or distilled water to a paste-like consistency.
3. Build dentin porcelain or opaqueous dentin porcelain around the cervical one-third. Condense slightly. If building a multiple-unit fixed-partial denture, apply the dentin porcelain mixture to the saddle area of the pontic and seat the framework on the model. Continue to apply the dentin porcelain until full contour has been achieved. Condense slightly.
4. Cutback the incisal one-third to allow room for

the enamel porcelains. If required, cutback to accentuate the lobe formation found in natural teeth. Be sure to remove adequate porcelain in the interproximal area.

5. Select the appropriate Natural Enamel Porcelain from **Table 1**. Apply the enamel porcelain to the incisal or occlusal one-third, overbuilding by approximately 10%.
6. Fire according to the recommended temperatures. A proper visual indicator is shiny with a very small amount of surface texture.

Advanced Porcelain Layering Procedures

1. Areas on the model that the porcelain will contact should be completely sealed with Ceramco PFZ Die Sealer. Other sealers or oils may leave a residue that will cause discoloration and should not be used.
2. Mix the porcelains with Ceramco PFZ Modeling Liquid, or distilled water to a paste-like consistency.
3. Build opaceous dentin/dentin modifier mixture with a 1:1 ratio around the cervical one third. Condense slightly. If building a multiple-unit fixed-partial denture, apply the opaceous dentin porcelain mixture to the saddle area of the pontic and seat the framework on the model.
4. Cutback the incisal one third to allow room for the enamel porcelains. If required, cutback to accentuate the lobe formation found in natural teeth. Be sure to remove adequate porcelain in the interproximal area.
5. If mamelons are required, place a small amount of mamelon porcelain on the dentin lobes. To simulate natural dentition, vary the length of the mamelons.
6. Select the appropriate Natural Enamel Porcelain from Table 1. Apply the enamel porcelain to the incisal or occlusal one-third, overbuilding by approximately 10%. If desired, layer the selected enamel porcelains in narrow vertical columns across the facial-incisal area to simulate the enamel rod formations found in natural teeth. Use Natural Enamel porcelains of different colors, or the clear enamel porcelain, or the Opal Enamel porcelains arranged in discrete parallel columns to create contrast. Condense lightly.
7. Remove the coping from the model. Apply a small amount of enamel porcelain to the mesial/distal contacts. Slightly condense the porcelain.

8. If building a multiple-unit fixed partial denture, use a thin bladed instrument and cut slightly into the interproximals. It is not necessary to cut all the way down to the liner.

9. Dry the build-up away from direct heat for at least five minutes before bringing it close to the open furnace entrance. Increase the drying time to ten minutes for large multiple-unit fixed partial dentures.
10. Fire according to the recommended temperatures. A proper visual indicator is shiny with a very small amount of surface texture.

Contouring & Cleaning

1. After firing, use diamond burs to refine the anatomy. Grinding too fast may cause excessive heat and crack the veneering porcelain. If necessary, use water to prevent excessive heat when grinding the Ceramco PFZ. If no porcelain additions are needed, the case may be glazed after thorough cleaning.
2. If additions are required or before glazing, lightly blast the porcelain surface with unrecycled 50-micron aluminum oxide at 20 psi. Clean the case with a steam cleaner or ultrasonic bath using distilled water.

Second Application of Dentin & Enamel Porcelains

Apply the additional porcelain using the same steps and techniques described for the initial application. Match the layering of the opal enamel porcelain over the dentin porcelain. Fire the porcelain using the firing parameters for the 2nd porcelain application.

Ceramco® PFZ Add-On Porcelains

The Ceramco PFZ Add-on porcelains are available in a Translucent enamel, and Tissue Tint shades Pink, Reddish-Pink, Salmon, and Dark. They fire at slightly lower than the regular dentin porcelain and are always fired under vacuum. The Translucent enamel porcelain may be used for low temperature additions, and the tissue colored porcelains would be used to simulate natural tissue.

1. Lightly blast the surface that is to be repaired with aluminum oxide at 20 psi. Then steam clean or clean in an ultrasonic cleaner in distilled water for 5 minutes.
2. Mix the porcelains with Ceramco PFZ Modeling Liquid, or distilled water to a paste-like consistency.

- 3. Apply the mixture to the desired area.
- 4. Fire according to the recommended firing temperatures.
- 3. Mix the stains with the stain liquid to a thin consistency and apply the mixture to the porcelain surface. Do not allow the mixture to puddle. Apply the stain until the desired effect is achieved.

Ceramco® PFZ Stains
Applying Ceramco PFZ Stains

- 1. Clean the porcelain surface by lightly blasting with 50 micron aluminum oxide at 20 psi, followed by cleaning with a steam cleaner or in an ultrasonic cleaner with distilled water.
- 2. Use the tip of the staining brush and spread a thin film of stain liquid and overglaze mixture over the surface to be stained. If you wish to use a natural glaze, use only the glaze medium. This is done to simulate a glaze and enable the ceramist to judge the necessary correction.
- 4. Fire according to recommended temperatures, using either Natural Glaze or Overglaze firing programs.
- 5. If stains are not required, the restoration can either be natural glazed, overglazed or mechanical polished using silicone wheels, porcelain polishing wheels, and flour pumice or aluminum oxide polishing paste. Diamond polishing paste may also be used but is not required.

Ceramco® PFZ Pressable Technique Table 1

Shade	Natural Enamel	Opal Effect Enamel	Margin Porcelain	Mamelon	Dentin Modifier	Shade
A1	Extra Light	Opal Light	A1	Yellow-Orange	A1	A1
A2	Light	Opal Light	A1	Yellow-Orange	A1	A2
A3	Light	Opal Light	A3	Red-Orange	A1	A3
A3.5	Light	Opal Medium	A3	Red-Orange	A3.5	A3.5
A4	Medium	Opal Medium	A4	Red-Orange	A3.5	A4
B1	Extra Light	Opal White	B2	Yellow-Orange	B1	B1
B2	Light	Opal Light	B2	Yellow-Orange	B1	B2
B3	Medium	Opal Light	A3	Red-Orange	B3	B3
B4	Medium	Opal Light	A4	Red-Orange	B3	B4
C1	Light	Opal Light	C1	Yellow-Orange	C1	C1
C2	Medium	Opal Light	C1	Yellow-Orange	C1	C2
C3	Medium	Opal Medium	C3	Red-Orange	C3	C3
C4	Medium	Opal Medium	C3	Red-Orange	C3	C4
D2	Extra Light	Opal Light	C1	Yellow-Orange	A1	D2
D3	Medium	Opal Light	A3	Yellow-Orange	B1	D3
D4	Medium	Opal Light	A3	Red-Orange	C3	D4

Ceramco® PFZ Press Firing Temperatures

FIRING PROGRAMS IN °C	Time (min)			Time (min)		Vacuum	Set Temp			Temperature		
	Dry	Pre Heat	Vac Hold	Hi Temp Hold	Cool	Set Point (in. Hg)	Idle	Hi Temp	Vac Start	Vac Stop	Heat Rate °C/min	Night
Margin/Margin Modifier	5	5	0	15 sec	0	29	450	970	450	970	60	100
Pressable Liner, 1st Bake	5	5	0	1 min	0	29	450	930	450	900	60	100
Pressable Liner, 2nd Bake	5	5	0	0	0	29	450	920	450	900	60	100
Op Dentin/Dentin/Mod/Enamels, 1st Bake	5	5	0	15 sec	0	29	450	900	450	900	60	100
Op Dentin/Dentin/Mod/Enamels, 2nd Bake	5	5	0	0	0	29	450	890	450	890	60	100
Natural Glaze/Stains	5	5	0	30 sec	0	0	450	870			60	100
Overglaze/Stains	5	5	0	30 sec	0	0	450	850			60	100
Add-On	5	5	0	0	0	29	450	860	450	860	60	100

FIRING PROGRAMS IN °F	Time (min)			Time (min)		Vacuum	Set Temp			Temperature		
	Dry	Pre Heat	Vac Hold	Hi Temp Hold	Cool	Set Point (in. Hg)	Idle	Hi Temp	Vac Start	Vac Stop	Heat Rate °F/min	Night
Margin/Margin Modifier	5	5	0	15 sec	0	29	842	1778	842	1778	108	212
Pressable Liner, 1st Bake	5	5	0	1 min	0	29	842	1706	842	1652	108	212
Pressable Liner, 2nd Bake	5	5	0	0	0	29	842	1688	842	1652	108	212
Op Dentin/Dentin/Mod/Enamels, 1st Bake	5	5	0	15 sec	0	29	842	1652	842	1652	108	212
Op Dentin/Dentin/Mod/Enamels, 2nd Bake	5	5	0	0	0	29	842	1634	842	1634	108	212
Natural Glaze/Stains	5	5	0	30 sec	0	0	842	1598			108	212
Overglaze/Stains	5	5	0	30 sec	0	0	842	1562			108	212
Add-On	5	5	0	0	0	29	842	1580	842	1580	108	212

Ceramco® PFZ Press Pressing Chart

Preheat the plungers only. Do not preheat the ingots.

PRESSING PROGRAMS IN °C	Low Temp °C	*High Temp °C	Heat Rate °C	High Hold Temp	Pressing Time 100g Ring Bench Cool	Pressing Time 200g Ring Bench Cool	Re-Press Time	Pressure	Vacuum
Finesse Press	700	940	60	20 min	10	20	na	4.25 bars	Full
Multimat Touch & Press	700	940	60	20 min	10	20	na	2.7 bars	Full
Empress EP 500	700	940	60	20 min	na	na	3	5.0 bars	Full
QEX CeramPress	700	940	60	20 min	10	20	na	4.25 bars	Full
Pro-Press 100	700	940	60	20 min	10	20	na	4.25 bars	Full

*Due to variance in furnaces, high temperatures may need to be adjusted slightly.

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