FINESSE®

LOW-FUSING
PORCELAIN

QUICK START INSTRUCTIONS

UNITED STATES

DENTSPLY Ceramco, Six Terri Lane, Burlington, NJ 08016, USA

Tel: 609-386-8900 • 800-487-0100

Fax: 609-386-8282

www.dentsply.com

EU-Rep. DENTSPLY Detrey GmbH,
Rodenbacher Chaussee 4, D-63457 Hanau-Wolfgang

DENTSPLY Canada
161 Vinyl Court, Woodbridge, Ont. L4L 4A3
905-851-6060

DENTSPLY

CERAMCO

PC 090172

CE 0086

ISO Certified 9001

FINESSE®

LOW-FUSING
PORCELAIN SYSTEM

© 2003 Ceramco Inc. All Rights Reserved. Printed in USA.

Rev. 1 03/03

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

FINESSE® QUICK START INSTRUCTIONS

PHOTOS

Figure 1.	First paste opaque application
Figure 2.	First paste opaque firing
Figure 3.	Second paste opaque application
Figure 4.	Second paste opaque firing
Figure 5.	First powder opaque application
Figure 6.	First powder opaque firing
Figure 7.	Second powder opaque application
Figure 8.	Second powder opaque firing
Figure 9.	Dentin build-up, full contour
Figure 10.	Dentin cut-back
Figure 11.	Opal enamel build-up, full contour
Figure 12.	Dentin fired
Figure 13.	Opaque application for Shoulder Porcelain
Figure 14.	First Shoulder Porcelain application
Figure 15.	Removing the Shoulder Porcelain from the die
Figure 16.	Second Shoulder Porcelain application
Figure 17.	Completed Porcelain Margin
Figure 18.	Overglaze application
Figure 19.	Polishing the porcelain surface
Figure 20.	Completed Finesse PFM

THE FINESSE SYSTEM

The Finesse porcelain system has been engineered to provide exceptional handling, thermal stability and esthetic properties. Finesse is a low-fusing porcelain system which polishes easily and is kinder against the opposing dentition than standard dental porcelain systems. The entire system, from the paste opaques through the dentins, opal enamels and surface stains, imparts the fluorescence of natural dentition. The true, life-like opalescence and fluorescence of the porcelains impart unsurpassed vitality.

The Finesse powder porcelains contain colorants for easy identification in porcelain application. The colorants disappear during firing.

INDICATIONS

The Finesse porcelains may be used for porcelain fused-to metal fixed prosthodontics using conventional high-fusing ceramic alloys.

CONTRAINDICATIONS

Only the indications listed above are suitable. These porcelains should not be used over aluminous porcelains or Dicor® glass-ceramic core materials because of thermal incompatibility.

PRECAUTIONS

Alloy Selection

Finesse porcelain is compatible with high noble, noble and predominantly base alloys. Alloys which generally fall into this category include high gold/platinum/palladium, high gold/platinum, 50% gold/ platinum/palladium, palladium, palladium/ low gold, palladium/ low gold/low silver and base alloys. Be sure to consult the alloy manufacturer for alloy compositions and coefficient of thermal expansion data to determine if the alloy is thermally compatible with Finesse porcelain.

Porcelain Furnace Programming

Finesse porcelains fire at significantly lower temperatures than high-fusing ceramic systems. Be sure to consult the firing schedule. Older porcelain furnaces may not be calibrated at the Finesse

porcelain firing temperatures. Fire a sample tab of Finesse porcelain before starting your first case. Be sure to obtain the correct visual indicator. If required, make necessary adjustments with all Finesse firing temperatures.

Sagger Trays

DO NOT USE SOLID SAGGER TRAYS.

During the firing cycle, solid ceramic sagger trays absorb more heat than honeycomb trays. Use of solid ceramic sagger trays may cause the Finesse porcelains to be underfired and appear chalky or porous. Use the honeycomb style trays with metal firing pins for the best results with the Finesse low-fusing system.

Liquids

USE FINESSE LIQUIDS ONLY.

Caution should be used in selecting modeling liquids. Finesse Liquids are recommended for their ability to improve the handling and carving of Finesse porcelain, to reduce slumping during build-up and completely burn-out at low temperatures. Finesse Liquids are specially formulated to completely burn-out at the low firing temperatures for Finesse porcelain. Liquids other than Finesse Liquids may not burn-out completely and may consequently make the porcelain gray.

Finesse Carving Liquid, Opaque Modifier Liquid, Margin Liquid and Stain Liquid should be used with the corresponding porcelains.

Brushes and Instruments

Use clean brushes and instruments. Rinse with clean water. Contamination with high-fusing porcelain will cause white spots in the fired porcelain.

STORAGE

Caution: THIS PORCELAIN IS SENSITIVE TO HUMIDITY.

Store porcelain in bottles with the lid tightly sealed. Do not leave bottles open. Avoid high heat, intense sunlight and high humidity. Shake well before use. The lot number is printed on each bottle; please refer to the lot number in all correspondence.

Mixed porcelain should not be stored for more than one days. The use of humidior-type trays for long-term storage is not recommended; white spots result from repeated mixing or wet storage of porcelain. If the tray or liquid should become contaminated, the porcelain may become gray and appear lower in value.

APPLICATION

Paste Opaque Application

1. Oxidize or degas the coping according to the alloy manufacturer's instructions.
2. Dispense 5-20 mm of paste onto a pallet. Using a brush, apply the opaque to the alloy in a thin, even layer. **(Figure 1)** A slight

vibration may be used to smooth out the opaque.

3. Dry in a warm area several minutes and then fire according to the schedule provided.
4. Check to see that the correct visual indicator has been achieved: moderate sheen. **(Figure 2)**
5. Apply a second coat of paste opaque porcelain to the alloy, being sure the metal is completely masked. If opaque modifiers are required, apply them at this time. **(Figure 3)**
6. Sprinkle the crystals on the second layer of the paste opaque only. Fire the restoration according to the schedule provided.
7. Check to see that the correct visual indicator has been achieved: sandpaper or emery board appearance. **(Figure 4)**

Powder Opaque Application

1. Oxidize or degas the alloy according to the alloy manufacturer's instructions.
2. Mix the powder opaque porcelain with Finesse Carving Liquid or distilled water to a thin paste-like consistency.
3. Paint the opaque on the alloy in a thin, even layer. Condense slightly. **(Figure 5)**
4. Dry in a warm area for several minutes and then fire according to the schedule provided.

5. Check to see that the correct visual indicator has been achieved: slight sheen. **(Figure 6)**
6. Paint a second coat of powder opaque on the alloy, being sure the metal is completely masked out. If opaque modifiers are required, apply them at this time. **(Figure 7)**
7. Condense slightly and fire according to the schedule provided.
8. Check to see that the correct visual indicator has been achieved: slight sheen. **(Figure 8)**

Basic Dentin/Opal Enamel Build-Up Technique

1. Seat the coping on the model.
2. Place the dentin, opaceous dentin, dentin modifier porcelain and an appropriate opal enamel porcelain (see page 6) on the palette. Mix with Finesse Carving Liquid or distilled water to a paste-like consistency.
3. Using a brush or an instrument, apply the dentin and/or opaceous dentin and dentin modifier porcelain until full contour has been achieved. Slightly condense the porcelain during the application. **(Figure 9)**
4. Remove the dentin porcelain from the incisal one-third. Slightly taper the “cutback” down to the gingival one-third. Be sure the dentin porcelain is covering all opaqued surfaces.

Note: Only a thin layer of opal

enamel porcelain is required. **(Figure 10)**

5. Apply the natural enamel or opal enamel porcelain to the incisal or occlusal one-third, overbuilding by approximately 10%. Condense slightly. **(Figure 11)**
6. Remove the coping from the model. Apply a small amount of opal enamel porcelain to the mesial/distal contacts. Slightly condense the porcelain.
7. Fire according to the schedule provided.
8. Check to see that the correct visual indicator has been achieved: shiny with a slightly textured appearance. **(Figure 12)**

Contouring & Cleaning

1. Use fine grit, non-contaminating stones, discs, or carbide and diamond burs to refine the anatomy. If no porcelain additions are needed, the case may be glazed or polished after thorough cleaning.
2. If additions are required or before glazing, lightly blast the porcelain surface with non-recycled 50 micron aluminum oxide at 40 psi. Clean the porcelain with a steam cleaner or ultrasonic bath using distilled water.

Second Application of Dentin and Opal Enamel Porcelains

1. Apply the additional porcelain using the same steps and techniques described for the initial

application. Match the layering of the opal enamel porcelain over the dentin porcelain.

2. Fire according to the schedule provided.

Tooth Margin Preparation

1. A 90° shoulder preparation or deep chamfer is recommended.
2. Normal tooth reduction (1.5mm) is recommended.

Opaque Porcelain Procedure

1. Apply and fire the Finesse paste opaque or Finesse powder opaque porcelain using your normal Finesse porcelain procedures. (**Figure 13**)

Die Sealer Liquid Application

1. Using the Finesse die sealer, apply a thin, even layer of sealer over the die margin area.
2. Blow off any excess die sealer and allow the liquid to thoroughly dry.

Die Release Liquid Application

1. Using the Finesse die release, apply a thin, even layer of die release over the margin area.
2. Blow off the excess die release.
3. Apply and allow to dry several times until the die has a very slight sheen.

First Margin Porcelain Application

1. Mix the Finesse margin liquid with the porcelain to a thick, creamy consistency. This mixture will air harden in approximately 5 minutes. If additional working time is required, mix in a small amount of distilled water.
2. Place the porcelain mixture in the gingival area, pushing the mixture down to, but not over the gingival margin. Do not overbuild. (**Figure 14**)
3. Allow the Finesse margin porcelain to dry completely (at least 5 minutes). The porcelain surface will be chalky and hard.
4. Carefully lift the coping from the die. (**Figure 15**)
Note: If inspection reveals the need for adjustment, carefully return the coping to the die and add fresh material. It is not necessary to remove the unfired porcelain, according to the alloy manufacturer's instructions.
5. Fire the restoration according to the margin porcelain firing schedule provided.

Second Margin Porcelain Application

1. Be sure the die is clean and reapply the die release. Allow to dry.
2. Make a new mix of margin porcelain and Finesse margin liquid. With the coping off the die,

- use a brush or instrument to apply a small amount of margin porcelain to the gingival margin.
- While the margin porcelain mixture is wet, firmly place the coping on the die and adjust the margin porcelain as necessary. (Figure 16)
 - Allow the margin porcelain to dry completely (at least 5 minutes). Remove the coping from the die.
 - Fire the restoration according to the margin firing schedule provided. (Figure 17)
- glaze or Finesse stains, mix with Finesse Stain Liquid to a thin consistency. Apply the mixture to the porcelain surface. Fire according to the schedule provided. (Figure 18)
- If desired, polish the porcelain after glazing. Porcelain polishing wheels, aluminum oxide polishing pastes, felt wheels or very fine-grit pumice compounds can be used for this technique. Diamond polishing paste may also be used but is not required to achieve a high luster. (Figure 19)

STAINING AND GLAZING

- Use a natural glaze cycle or apply Finesse Glaze.

- Completed case. (Figure 20)

SUGGESTED OPAL ENAMELS

Shade	Primary Option	Secondary Option	Enamel Effect Color
A0/i1	Natural White	Opal White	Opal Pink
A1	Natural Light	Opal Light	Opal Pink
A2	Natural Light	Opal Light	Opal Pink or Opal Yellow
A3	Natural Light	Opal Light	Opal Orange or Opal Pink
A3.5	Natural Medium	Opal Medium	Opal Orange
A4	Natural Medium	Opal Medium	Opal Orange
B0/i2	Natural White	Opal White	Opal Yellow
B1	Natural White	Opal White	Opal Yellow
B2	Natural Light	Opal Light	Opal Yellow
B3	Natural Light	Opal Light	Opal Orange
B4	Natural Light	Opal Light	Opal Orange
C1	Natural Light	Opal Light	Opal Yellow
C2	Natural Light	Opal Light	Opal Orange
C3	Natural Medium	Opal Medium	Opal Orange
C4	Natural Medium	Opal Medium	Opal Orange
D2	Natural Light	Opal Light	Opal Pink
D3	Natural Light	Opal Light	Opal Pink
D4	Natural Light	Opal Light	Opal Orange

PRODUCT AVAILABILITY

Finesse is a complete system of low-fusing porcelains. Finesse porcelain was developed to allow ease of use in creating esthetic restorations.

The complete Finesse system includes:

Powder Opaque Porcelains	18 Shades
Powder Opaque Modifiers	9 Shades
Paste Opaque Porcelains	18 Shades + 1 Crystal
Paste Opaque Modifiers	9 Shades
Opaque Correction Pastes	3 Shades
Opaceous Dentin Porcelains	18 Shades
Dentin Porcelains	18 Shades
Dentin Modifiers	18 Shades + 8 Primary Colors
Dentin Correction Porcelains	3 Shades
Enamel Porcelains	15 Shades
Natural Enamels	4 Shades
Mamelon Porcelains	3 Shades
Tissue Tint Porcelains	2 Shades
Margin Porcelains	18 Shades
Margin Modifiers	6 Shades
Stains	18 Shades
Shade Fans	9 Fans
Glaze	
Stain Palette	
Opaque Mixing Palette	
Syringe Stand	
Paste Opaque Brush	
Crystal Collection Bowl	
Opaque Modifier Liquid	
Margin Liquid	
Carving Liquid	
Stain Liquid	
Die Release	
Die Sealer	

NOTE: To achieve optimum results, observe the guidelines and instructions provided in this brochure. Separate instruction manuals are available for the components of this system.

FINESSE® FIRING REFERENCE CHARTS

FIRING SCHEDULES FOR FINESSE PORCELAIN (°C) Use honeycomb sagger trays with metal pins.

Firing Program in °C	Time (min.)			Time (min.)			Vacuum		Set Temp.		Temperature		
	Dry	Pre-Heat	Vac. Hold	Hi Temp Hold	Hold	Cool	Set Point (in. Hg)	Idle	Hi Temp	Vac. Start	Vac. Stop	Heat Rate (°C/min.)	Night
Powder Opaque	3	3	0	1	0	0	29	450	800	450	800	90	95
Paste Opaque	5	3	0	0	0	0	29	450	790	450	790	90	95
Opaque Correction Opaqueous Dentin, Dentin Modifiers Dentin Tissue Tint	5	5	0	0.5	0	0	29	450	760	450	760	35	95
1st Firing Dentin/ Enamel	5	5	0	0.5	0	0	29	450	760	450	760	35	95
2nd Firing Dentin/Enamel	5	5	0	0	0	0	29	450	750	450	750	35	95
Dentin Correction	5	5	0	0	0	0	29	450	730	450	710	55	95
Margin Margin Modifiers	3	7	0	0	0	0	29	675	770	675	770	35	95
Natural Glaze	3	3	0	0	0	0	0	450	750	na	na	70	95
Glaze/Stain	3	3	0	0	0	0	0	450	750	na	na	70	95

Large cases and/or multiple units may require a higher temperature or a lower heat rate. Do not fire porcelain with carbon blocks. If the results achieved with these schedules do not conform to your needs, refer to the Trouble Shooting Section.

FIRING SCHEDULES FOR FINESSE PORCELAIN (°F) Use honeycomb sagger trays with metal pins.

Firing Program in °F	Time (min.)			Time (min.)			Vacuum		Set Temp.		Temperature		
	Dry	Pre-Heat	Vac. Hold	Hi Temp Hold	Hold	Cool	Set Point (in. Hg)	Idle	Hi Temp	Vac. Start	Vac. Stop	Heat Rate (°F/min.)	Night
Powder Opaque	3	3	0	1	0	0	29	842	1472	842	1472	162	200
Paste Opaque	5	3	0	0	0	0	29	842	1454	842	1454	162	200
Opaque Correction Opaqueous Dentin, Dentin Modifiers Dentin Tissue Tint	5	5	0	0.5	0	0	29	842	1400	842	1400	63	200
1st Firing Dentin/ Enamel	5	5	0	0.5	0	0	29	842	1400	842	1400	63	200
2nd Firing Dentin/Enamel	5	5	0	0	0	0	29	842	1382	842	1382	63	200
Dentin Correction	5	5	0	0	0	0	29	842	1346	842	1310	99	200
Margin Margin Modifiers	3	7	0	0	0	0	29	1247	1418	1247	1418	63	200
Natural Glaze	3	3	0	0	0	0	0	842	1382	na	na	126	200
Glaze/Stain	3	3	0	0	0	0	0	842	1382	na	na	126	200

Large cases and/or multiple units may require a higher temperature or a lower heat rate. Do not fire porcelain with carbon blocks. If the results achieved with these schedules do not conform to your needs, refer to the Trouble Shooting Section.

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