

Perfection in Porcelain.

DUCERAGOLD



Duceragold Duceram Plus Excelsior LFC AllCeram Finesse Ceramco 3

Duceragold — the premium porcelain that's compatible with a variety of alloys.

- Advanced formulation for the ultimate in esthetics and function
- Assures a dense, smooth restoration surface which resists plaque accumulation
- Forms a self-repairing protective hydroxyl layer in the oral environment
- Provides high oral acid abrasion resistance virtually equal to natural dentition
- The original wear-friendly porcelain with the beauty of Light Dynamics
- Compatible with type IV gold alloys such as Degunorm[®], Degudent[®] LS, Stabilor[®] LS, or Sonoma[®] and Cergo[®] All-Ceramic
- Low fusing Ceramic with high firing temperature of 800°C (Dentins)
- CTE Range: 16.2 - 17.1

For more information about Duceragold Porcelain, call DENTSPLY Ceramco at 800-487-0100 or visit our website at www.dentsply.com

Perfection in Porcelain.

DUCERAGOLD



Duceragold Assortments

THE PROFESSIONAL SET – PC 60090071

| | |
|----------|-----------------------|
| 4 x 20g | Opal Incisals 57 – 60 |
| 5 x 20g | Opal Incisals 10 – 50 |
| 4 x 20g | Flu Dentin Porcelains |
| 10 x 20g | Modifiers |
| 1 x 50ml | SD Modeling Liquid |
| 1 | Shade Indicator |

THE DENTIN PORCELAIN SET – PC 60090041

| | |
|----------|----------------------|
| 16 x 20g | Dentins A1 – D4 |
| 3 x 20g | 3D Dentin Porcelains |
| 1 x 50ml | SD Modeling Liquid |
| 1 | Shade Indicator |

THE OPACEOUS DENTIN PORCELAIN SET – PC 60090031

| | |
|----------|---------------------------|
| 16 x 20g | Opaceous Dentins A1 – D 4 |
| 1 x 50ml | SD Modeling Liquid |
| 1 | Shade Indicator |

BLEACH SET – PC 60090091

| | |
|----------|-----------------------------|
| 3 x 20g | Bleach Dentins 1 – 3 |
| 1 x 2ml | Bleach Paste Opaque |
| 1 x 20g | Transparent Porcelain |
| 1 x 20g | Bleach Shoulder Margin |
| 1 x 15ml | SD Liquid |
| 1 x 15ml | SM Liquid |
| 1 x 10ml | DUCERA SEP Separating Agent |
| 1 | Paste Opaque Brush |
| 1 | Shade Indicator |

THE INCISAL PORCELAIN SET – PC 60090051

| | |
|----------|---|
| 4 x 20g | Incisals 57 – 60 |
| 3 x 20g | Transparent Porcelains (TC, T, TO) |
| 3 x 20g | Creative Incisal Porcelains (reddish, blueish, yellowish) |
| 2 x 20g | Gingival Porcelains (light, dark) |
| 1 x 20g | Add-on Porcelain |
| 1 x 20g | Glaze |
| 1 x 50ml | SD Modeling Liquid |
| 1 x 10ml | DUCERA SEP Separating Agent |
| 1 x 15ml | Stain Liquid |
| 1 | Shade Indicator |

THE LIGHT DYNAMICS ASSORTMENT – PC 60090061

| | |
|----------|--------------------|
| 11 x 20g | True Opal Incisals |
| 4 x 20g | Twin Dentins |
| 1 x 50ml | SD Modeling Liquid |
| 4 x 20g | Bright Body Flu |
| 1 | Shade Indicator |

SHOULDER MARGIN SET – PC 60090081

| | |
|----------|--|
| 8 x 20g | Shoulder Margin Material 1 – 7 and Flu Fluorescent |
| 1 x 50ml | SD Quick Liquid |
| 1 | Shade Indicator |

POWDER OPAQUE SET – PC 60090011

| | |
|----------|------------------|
| 1 x 20g | Bond |
| 16 x 20g | Opagues A1 – D 4 |
| 7 x 20g | Opaque Modifiers |
| 1 x 50ml | OL Liquid |
| 1 | Shade Indicator |

PASTE OPAQUE SET – PC 60090021

| | |
|----------|------------------------|
| 16 x 2ml | Paste Opaque A1 – D4 |
| 7 x 2ml | Paste Opaque Modifiers |
| 1 x 2ml | Redox Multi |
| 2 | Paste Opaque Brushes |
| 1 | Shade Indicator |

DUCERAGOLD TEST SET WITH CERGO PELLETS – PC 60090100

| | |
|----------|---------------------------------------|
| 1 x 4g | Add-on Porcelain |
| 1 x 4g | 3D Dentin Porcelain |
| 2 x 4g | Dentin Porcelains A2 & A3 |
| 2 x 4g | Bright Body Flu |
| 2 x 4g | Twin Dentins |
| 6 x 4g | True Opal Incisals |
| 2 x 1ml | Paste Opagues A2 & A3 |
| 1 x 1ml | Redox Multi |
| 1 x 15ml | SD Liquid |
| 1 x 15ml | SD Form Liquid |
| 1 | Paste Opaque Brush |
| 1 Each | Cergo Pellets A2, A3, A3.5, S1, S3, T |
| 2 Each | Cergo Pellets S2 & TR |

DENSPLY
CERAMCO

The chemistry of a great smile.

Six Terri Lane Burlington, NJ 08016
1-800-487-0100 www.dentsply.com